








Cushions for pressure relief of body and arms



all up[®]


Name		Size in cm.	Item no.	Application	About the device
All up - Multi Small		50 x 25 x 8	70003	Positioning cushion with 2 chambers. Multifunctional. Can be used for both supine and sitting positions. (between legs, underneath the back, providing excellent support under arms and between the knees.).	
All up - Multi Small Silent & Soft		50 x 25 x 8	70105	Positioning cushion with 2 chambers. Multifunctional. Can be used for both supine and sitting positions. (between legs, underneath the back, providing excellent support under arms and between the knees.), Extra soft and silent pillow.	
All up - Multi Medium		50 x 37 x 8	70004	Positioning cushion with 3 chambers. Multifunctional. Especially useful for a supine position (between legs, back), and possibly as support under the feet. Extra soft and silent pillow.	
All up - Multi Large		50 x 50 x 8	70005	Positioning cushion with 4 chambers. Multifunctional. Can be folded/rolled up. For use under the leg or arm, comes with an extra Velcro strap.	
All up - Multi Max		50 x 50 x 20	70012	Positioning cushion, particularly useful for hip fracture patients and when the lower limb needs to be elevated above heart level.	
All up - Universal		50 x 50 x 8	70010	Positioning cushion, typically for the back, sacrum and tailbone in adjustable bed and the lateral positions between the legs. Can be folded if a higher cushion is desired.	
All up - Universal Long		55 x 70 x 8	70015	Positioning cushion, typically for the lateral positions between the legs, longer model which protects the knees and malleolus. Can be folded in half if a higher cushion is desired.	
All up - Universal Slim		45 x 75 x 8	70016	Positioning cushion, typically for lateral positions between the legs, longer and narrower model which ensures a distance between the knees in the lateral position. Can be folded if a higher cushion is desired.	
All up - Head		28 x 30 x 5	70011	Cushion, particularly useful for the preventing and treating of pressure sores on the ears and neck, and is also useful for coma patients.	
All up - Head Silent & Soft		28 x 30 x 5	70074	Cushion, particularly useful for the preventing and treating of pressure sores on the ears and neck, and is also useful for coma patients. Extra soft and silent pillow.	
All up - Wing		50 x 85 x 8	70006	Relief cushion for lower limbs, used for bed-bound patients. The cushion is easily fixed to the mattress.	
All up - Arm Rest		55 x 25 x 8	70027	Positioning cushion for supporting the arms.	
All up - Elbow		30 x 20 x 10	70028	Positioning cushion for supporting the elbow.	
All Up - Elbow Medium		23 x 20 x 10	70048	Positioning cushion for supporting the elbow.	



Name		Size in cm.	Item no.	Application	About the device
All up - Lap		80 x 60 x 10	70025	Positioning cushion, crescent. Advantageous for wheelchair users for supporting the arms.	
All up - Lap Max		80 x 60 x 20	70026	Positioning cushion, crescent. With extra height. Advantageous for wheelchair users for supporting the arms.	
All Up - Elbow & Heel		16 x 19 x 10	70049	Cushion for relieving heel or elbow.	

all^{up} kids^{care}

All Up- Kids Care Head		16 x 18 x 5	70040	Small pillow. Particularly useful for children under 1 year being ventilated	
Kids Care Head Silent & Soft		16 x 18 x 5	70075	Small pillow. Particularly useful for children under 1 year being ventilated	



Cushions for special positioning




Name		Size in cm.	Item no.	Application	About the device
All up - Belly		65 x 75 x 15	70023	Special positioning cushions with particular relevance for radiology departments to elevate the chest during examinations or ventilating in prone position	
All up - Kidney		50 x 20 x 20	70014	Special positioning cushion with particular relevance for radiology departments for supporting the kidney/pelvic area.	
All up - Wedge 30°		50 x 40 x 15	70024	Positioning cushion for the lateral position. Ensures a stable position at 30°, and the wedge shape is fixed by a flap that is pushed underneath the patient.	
All up - Multi Wedge		50 x 23 x 23	70053	Slanted cushion that elevates the patient's arm, which is stable in the pillow. Can be used for lymphedema patients (venous catheter in hand)	
All up - Head Intubate DUO White		32 x 30 x 18	70103	Cushion with tube openings for ventilating in the prone position.	
All Up - Foot sole		86 x 37 x 5	70080	Cushion for relieving the foot soles. Mounted on the bed rail at the foot end of the bed. Fits all beds.	

PRONE kit

Prone Kit incl. manual pump	 	The set contains, item no.: 70023 70103	70115	Set consisting of two cushions: All Up Head intubate DUO White and All Up Belly is designed to be used in patients who are intubated in the prone position. All Up Belly can be placed on the patient before the patient is turned over and is pumped up when the patient is turned back in the prone position. The set contains 3 pcs. of each cushion. The Cushions can be used individually.	
-----------------------------	------	---	-------	---	--

COVID19 relief kit

Covid 19 Relief Kit incl. manual pump	 	The set contains, item no.: 70023 70103 70011 70010 70003 70007 70002	70104	Set consisting of 7 cushions: All Up Head Intubate DUO White, All Up Belly, All Up Head, All Up Universal, All Up Multi Small, Heel Up Medium and Heel Up Long. The set covers most positioning situations that may arise. Single-patient use ensures optimal hygiene.
---------------------------------------	--------------------------	--	-------	---





PRONE kit

Set consisting of two cushions:
All Up Head intubate DUO White
and All Up Belly is designed to be
used in patients who are intubated
in the prone position.

All Up Belly can be placed on the
patient before the patient is turned
over and is pumped up when the
patient is turned back in the prone
position.

The set contains 3 pcs.
of each cushion.

The cushions can be used
individually.




Seat All® has some important features that provide a number of advantages

**Drainage
of heat**

**Drainage
of moisture**

Lightweight



The outer surface is made of
nonwoven felt which provides
superior drainage of moisture and
heat from the skin, the felt also
ensures that the skin does not stick
to the surface of the cushion.

Seat All® can be wiped with a damp
cloth, alcohol or other suitable
disinfectant.